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
The screenshot displays the Shimadzu software interface for PHD analysis. The main window is titled 'PHD测定' (PHD Measurement). On the left, there are settings for '测定条件' (Measurement Conditions) and 'PHD测定条件' (PHD Measurement Conditions). The '测定条件' section includes '保存位置' (Save Location), '透光片' (Transmittance Plate), '检测器高压(V)' (Detector High Voltage), 'PHA(Low)', 'PHA(High)', '电压(V)' (Voltage), '电流(mA)' (Current), and 'PHD测定条件' (PHD Measurement Conditions) such as '启动级' (Start Level), '结束级' (End Level), '步速' (Step Speed), and '精扫时间(秒)' (Fine Scan Time). The 'PHD测定条件' section includes '启动级' (Start Level), '结束级' (End Level), '步速' (Step Speed), and '精扫时间(秒)' (Fine Scan Time).

A 'Group Select' dialog box is open, showing a list of samples with columns for '组名' (Group Name), '日期' (Date), and '备注' (Remarks). The '日期' column is highlighted, and a circled '3' is next to the '日期' header.

The 'PHD测定结果' (PHD Measurement Results) table is visible, showing columns for '元素' (Element), '分析线' (Analysis Line), '晶体' (Crystal), '检测器' (Detector), '检测器高压(V)' (Detector High Voltage), 'PHA(Low)', 'PHA(High)', '电压(V)' (Voltage), '电流(mA)', and '透光片' (Transmittance Plate). A circled '4' is next to the '检测器' column.

A graph shows the resulting spectrum with a peak labeled '5'.

The bottom part of the screenshot shows the 'PHD测定' (PHD Measurement) panel with a circled '1' next to the 'PHD测定' button. The 'PHD测定结果' (PHD Measurement Results) table is also visible, showing a circled '6' next to the 'PHD测定结果' button.

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