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
- [1] European Medicines Agency, Guideline on Bioanalytical Method Validation, 2011.
- [2] The International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use. Bioanalytical method validation and study sample analysis M10[S]. 2022.
- [3] 国家药典委员会. 中华人民共和国药典（四部）[S]. 北京: 中国医药科技出版社, 2025.
- [4] U.S. Department of Health and Human Services Food and Drug Administration. Bioanalytical Method Validation for Biomarkers Guidance for Industry[S]. 2025.
- [5] 叶潇, 郑爽, 林超, 等. 生物分析方法验证及生物样品分析各指导原则间的比较[J]. 中国药物评价, 2023, 40(06): 470-479.

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**附录 ICH M10 2020_M10 Bioanalytical
Method Validation and Study Sample
Analysis (节选)**

**M10 BIOANALYTICAL
METHOD VALIDATION
AND STUDY SAMPLE
ANALYSIS**
Guidance for Industry

**U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)
Center for Biologics Evaluation and Research (CBER)**

**November 2022
ICH**

M10 BIOANALYTICAL METHOD VALIDATION AND STUDY SAMPLE ANALYSIS

Guidance for Industry

*Additional copies are available from:
Office of Communications, Division of Drug Information
Center for Drug Evaluation and Research
Food and Drug Administration
10001 New Hampshire Ave., Hillandale Bldg., 4th Floor
Silver Spring, MD 20993-0002
Phone: 855-543-3784 or 301-796-3400; Fax: 301-431-6353
Email: druginfo@fda.hhs.gov*

*<https://www.fda.gov/drugs/guidance-compliance-regulatory-information/guidances-drugs>
and/or*

*Office of Communication, Outreach and Development
Center for Biologics Evaluation and Research
Food and Drug Administration
10903 New Hampshire Ave., Bldg. 71, Room 3128
Silver Spring, MD 20993-0002
Phone: 800-835-4709 or 240-402-8010
Email: ocod@fda.hhs.gov*

<https://www.fda.gov/vaccines-blood-biologics/guidance-compliance-regulatory-information-biologics/biologics-guidances>

**U.S. Department of Health and Human Services
Food and Drug Administration
Center for Drug Evaluation and Research (CDER)
Center for Biologics Evaluation and Research (CBER)**

**November 2022
ICH**

Contains Nonbinding Recommendations

FOREWORD

The International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use (ICH) has the mission of achieving greater regulatory harmonization worldwide to ensure that safe, effective, and high-quality medicines are developed, registered, and maintained in the most resource-efficient manner. By harmonizing the regulatory expectations in regions around the world, ICH guidelines have substantially reduced duplicative clinical studies, prevented unnecessary animal studies, standardized safety reporting and marketing application submissions, and contributed to many other improvements in the quality of global drug development and manufacturing and the products available to patients.

ICH is a consensus-driven process that involves technical experts from regulatory authorities and industry parties in detailed technical and science-based harmonization work that results in the development of ICH guidelines. The commitment to consistent adoption of these consensus-based guidelines by regulators around the globe is critical to realizing the benefits of safe, effective, and high-quality medicines for patients as well as for industry. As a Founding Regulatory Member of ICH, the Food and Drug Administration (FDA) plays a major role in the development of each of the ICH guidelines, which FDA then adopts and issues as guidance to industry.

Contains Nonbinding Recommendations

TABLE OF CONTENTS

I.	INTRODUCTION (1)	1
A.	Objective (1.1)	1
B.	Background (1.2)	1
C.	Scope (1.3)	2
II.	GENERAL PRINCIPLES (2)	3
A.	Method Development (2.1)	3
B.	Method Validation (2.2)	3
1.	<i>Full Validation (2.2.1)</i>	3
2.	<i>Partial Validation (2.2.2)</i>	4
3.	<i>Cross Validation (2.2.3)</i>	4
III.	CHROMATOGRAPHY (3)	4
A.	Reference Standards (3.1)	4
B.	Validation (3.2)	5
1.	<i>Selectivity (3.2.1)</i>	5
2.	<i>Specificity (3.2.2)</i>	6
3.	<i>Matrix Effect (3.2.3)</i>	7
4.	<i>Calibration Curve and Range (3.2.4)</i>	7
5.	<i>Accuracy and Precision (3.2.5)</i>	8
6.	<i>Carryover (3.2.6)</i>	9
7.	<i>Dilution Integrity (3.2.7)</i>	9
8.	<i>Stability (3.2.8)</i>	10
9.	<i>Reinjection Reproducibility (3.2.9)</i>	12
C.	Study Sample Analysis (3.3)	13
1.	<i>Analytical Run (3.3.1)</i>	13
2.	<i>Acceptance Criteria for an Analytical Run (3.3.2)</i>	14
3.	<i>Calibration Range (3.3.3)</i>	15
4.	<i>Reanalysis of Study Samples (3.3.4)</i>	15
5.	<i>Reinjection of Study Samples (3.3.5)</i>	16
6.	<i>Integration of Chromatograms (3.3.6)</i>	17
IV.	LIGAND BINDING ASSAYS (4)	17
A.	Key Reagents (4.1)	17
1.	<i>Reference Standard (4.1.1)</i>	17
2.	<i>Critical Reagents (4.1.2)</i>	17
B.	Validation (4.2)	18
1.	<i>Specificity (4.2.1)</i>	18
2.	<i>Selectivity (4.2.2)</i>	19
3.	<i>Calibration Curve and Range (4.2.3)</i>	19
4.	<i>Accuracy and Precision (4.2.4)</i>	20
5.	<i>Carryover (4.2.5)</i>	21
6.	<i>Dilution Linearity and Hook Effect (4.2.6)</i>	21
7.	<i>Stability (4.2.7)</i>	22

Contains Nonbinding Recommendations

C.	Study Sample Analysis (4.3)	23
1.	<i>Analytical Run (4.3.1)</i>	23
2.	<i>Acceptance Criteria for an Analytical Run (4.3.2)</i>	23
3.	<i>Calibration Range (4.3.3)</i>	24
4.	<i>Reanalysis of Study Samples (4.3.4)</i>	25
V.	INCURRED SAMPLE REANALYSIS (5)	26
VI.	PARTIAL AND CROSS VALIDATION (6)	27
A.	Partial Validation (6.1)	27
B.	Cross Validation (6.2)	28
VII.	ADDITIONAL CONSIDERATIONS (7)	29
A.	Methods for Analytes That Are Also Endogenous Molecules (7.1)	29
1.	<i>Quality Control Samples for Methods for Analytes That Are Also Endogenous Molecules (7.1.1)</i> ..	31
2.	<i>Selectivity, Recovery, and Matrix Effects for Methods for Analytes That Are Also Endogenous Molecules (7.1.2)</i>	31
3.	<i>Parallelism for Methods for Analytes That Are Also Endogenous Molecules (7.1.3)</i>	32
4.	<i>Accuracy and Precision for Methods for Analytes That Are Also Endogenous Molecules (7.1.4)</i>	32
5.	<i>Stability for Methods for Analytes That Are Also Endogenous Molecules (7.1.5)</i>	33
B.	Parallelism (7.2)	33
C.	Recovery (7.3)	33
D.	Minimum Required Dilution (7.4)	34
E.	Commercial and Diagnostic Kits (7.5)	34
F.	New or Alternative Technologies (7.6)	35
1.	<i>Dried Matrix Methods (7.6.1)</i>	35
VIII.	DOCUMENTATION (8)	36
A.	Summary Information (8.1)	36
B.	Documentation for Validation and Bioanalytical Reports (8.2)	37
	GLOSSARY	45

M10 BIOANALYTICAL METHOD VALIDATION AND STUDY SAMPLE ANALYSIS Guidance for Industry¹

This guidance represents the current thinking of the Food and Drug Administration (FDA or Agency) on this topic. It does not establish any rights for any person and is not binding on FDA or the public. You can use an alternative approach if it satisfies the requirements of the applicable statutes and regulations. To discuss an alternative approach, contact the FDA office responsible for this guidance as listed on the title page.

I. INTRODUCTION (1)²

A. Objective (1.1)

This guidance is intended to provide recommendations for the validation of bioanalytical methods for chemical and biological drug quantification and their application in the analysis of study samples. Adherence to the principles presented in this guidance will ensure the quality and consistency of the bioanalytical data in support of the development and market approval of both chemical and biological drugs.

The objective of the validation of a bioanalytical method is to demonstrate that it is suitable for its intended purpose. You may use an alternative approach from the recommendations in this guidance if it satisfies the requirements of applicable statutes and regulations. Alternative approaches should be supported by appropriate scientific justification. Applicants are encouraged to consult the regulatory authority(ies) regarding significant changes in method validation approaches when an alternate approach is proposed or taken. This guidance replaces the ICH draft guidance for industry *M10 Bioanalytical Method Validation*, issued on June 27, 2019.³

B. Background (1.2)

Concentration measurements of chemical and biological drug(s) and their metabolite(s) in biological matrices are an important aspect of drug development. The results of studies employing such methods contribute to regulatory decisions regarding the safety and efficacy of

¹ This guidance was developed within an Expert Working Group of the International Council for Harmonisation of Technical Requirements for Registration of Pharmaceuticals for Human Use (ICH) and has been subject to consultation by the regulatory parties, in accordance with the ICH process. This document was endorsed by the Regulatory Members of the ICH Assembly at *Step 4* of the ICH process on May 24, 2022.

² The numbers in parentheses reflect the organizational breakdown of the document endorsed by the Regulatory Members of the ICH Assembly at *Step 4* of the ICH process, May 24, 2022.

³ We update guidances periodically. For the most recent version of a guidance, check the FDA guidance web page at <https://www.fda.gov/regulatory-information/search-fda-guidance-documents>.

Contains Nonbinding Recommendations

drug products. It is therefore critical that the bioanalytical methods used are well characterized, appropriately validated, and documented in order to ensure reliable data to support regulatory decisions.

This guidance intends to facilitate development of drugs in accordance with the principles of the 3Rs (Reduce, Refine, Replace) for animal studies, where valid.

C. Scope (1.3)

This guidance describes the validation of bioanalytical methods and study sample analysis that are expected to support regulatory decisions. The guidance is applicable to the bioanalytical methods used to measure concentrations of chemical and biological drug(s) and their metabolite(s) in biological samples (e.g., blood, plasma, serum, other body fluids or tissues) obtained in nonclinical toxicokinetic (TK) studies conducted according to the principles of Good Laboratory Practice (GLP), nonclinical pharmacokinetic (PK) studies conducted as surrogates for clinical studies, and all phases of clinical trials, including comparative bioavailability/bioequivalence (BA/BE) studies, in regulatory submissions. Full method validation is recommended for the primary matrix intended to support regulatory submissions. Additional matrices should be validated as necessary.

For studies that are not submitted for regulatory approval or not considered for regulatory decisions regarding safety, efficacy, or labeling (e.g., exploratory investigations), applicants may decide on the level of qualification that supports their own internal decision-making.

The information in this guidance applies to the quantitative analysis by ligand binding assays (LBAs) and chromatographic methods such as liquid chromatography (LC) or gas chromatography (GC), which are typically used in combination with mass spectrometry (MS) detection.

For studies that are subject to GLP or Good Clinical Practice (GCP) requirements, the bioanalysis of study samples must conform to those requirements. The bioanalysis of biomarkers and bioanalytical methods used for the assessment of immunogenicity are not within the scope of this guidance.

In general, FDA's guidance documents do not establish legally enforceable responsibilities. Instead, guidances describe the Agency's current thinking on a topic and should be viewed only as recommendations, unless specific regulatory or statutory requirements are cited. The use of the word should in Agency guidances means that something is suggested or recommended, but not required.

Contains Nonbinding Recommendations

II. GENERAL PRINCIPLES (2)

A. Method Development (2.1)

The purpose of bioanalytical method development is to define the design, operating conditions, limitations, and suitability of the method for its intended purpose and to ensure that the method is ready for validation.

Before or during the development of a bioanalytical method, the applicant is encouraged to, if feasible, understand the analyte of interest (e.g., the physicochemical properties of the drug, in vitro and in vivo metabolism, preferential distribution between red blood cells and plasma, and protein binding) and consider aspects of any prior analytical methods that may be applicable.

Method development involves identifying the procedures and conditions involved with quantifying the analyte. Method development can include the characterization of the following bioanalytical elements: reference standards, critical reagents, calibration curve, quality control samples (QCs), selectivity and specificity, sensitivity, accuracy, precision, recovery, stability of the analyte, and minimum required dilution (MRD).

Bioanalytical method development does not require extensive record keeping or notation. Once the method has been developed, bioanalytical method validation proves that the method is suited to the analysis of the study samples.

If a problem is encountered with the method during the analysis of nonclinical or clinical study samples that requires that the analysis be stopped, any changes to the method and the rationale should be documented.

B. Method Validation (2.2)

1. Full Validation (2.2.1)

Bioanalytical method validation is important to ensure the acceptability of assay performance and the reliability of analytical results. A bioanalytical method is defined as a set of procedures used for measuring analyte concentrations in biological samples. A full validation of a bioanalytical method should be performed when establishing a bioanalytical method for the quantification of an analyte in clinical and in applicable nonclinical studies. Full validation should also be performed when implementing an analytical method that is reported in the literature and when a commercial kit is repurposed for bioanalytical use in drug development. Usually, one analyte has to be determined, but on occasion it may be appropriate to measure more than one analyte. This may involve two different drugs, a parent drug with its metabolites or the enantiomers or isomers of a drug. In these cases, the principles of validation and analysis apply to all analytes of interest. For chromatographic methods, a full validation should include the following elements, unless otherwise justified: selectivity, specificity, matrix effect, calibration curve (response function), range (lower limit of quantification (LLOQ) to upper limit

Contains Nonbinding Recommendations

of quantification (ULOQ)), accuracy, precision, carryover, dilution integrity, stability, and reinjection reproducibility.

For LBAs, the following elements should be evaluated, unless otherwise justified: specificity, selectivity, calibration curve (response function), range (LLOQ to ULOQ), accuracy, precision, carryover, dilution linearity, and stability. If necessary, parallelism can be conducted when appropriate study samples are available.

The assessments that are performed during validation should be relevant to the sample analysis workflow. The matrix used for bioanalytical method validation should be the same as the matrix of the study samples, including anticoagulants and additives. In cases in which it may be difficult to obtain an identical matrix to that of the study samples (e.g., rare matrices such as tissue, cerebrospinal fluid, bile or in cases where free drug is measured), surrogate matrices may be acceptable for analytical method validation.

The choice of surrogate matrix should be scientifically justified. Matrix differences within species (e.g., age, ethnicity, sex) are generally not considered different when validating a method.

A specific, detailed, written description of the bioanalytical method and validation procedure should be established *a priori*. This description may be in the form of a protocol, study plan, report, notebook, or Standard Operating Procedure (SOP).

2. *Partial Validation (2.2.2)*

Modifications to a fully validated analytical method may be evaluated by partial validation. Partial validation can range from as little as one accuracy and precision determination to a nearly full validation (refer to section VI.A (6.1)). The items in a partial validation should be determined according to the extent and nature of the changes made to the method.

3. *Cross Validation (2.2.3)*

Cross validation is required to demonstrate how the reported data are related when multiple bioanalytical methods and/or multiple bioanalytical laboratories are involved (refer to section VI.B (6.2)).

III. CHROMATOGRAPHY (3)

A. Reference Standards (3.1)

During method validation and the analysis of study samples, a blank biological matrix is spiked with the analyte(s) of interest using solutions of reference standard(s) to prepare calibration standards and QCs. Calibration standards and QCs should be prepared from separate stock solutions. However, calibration standards and QCs may be prepared from the same stock

Contains Nonbinding Recommendations

solution provided the accurate preparation and stability of the stock solution should have been verified.

A suitable internal standard (IS) should be added to all calibration standards, QCs, and study samples during sample processing. The absence of an IS should be justified.

It is important that the reference standard is well characterized and the quality (e.g., purity, identity) of the reference standard and the suitability of the IS is ensured, as the quality will affect the outcome of the analysis and, therefore, the study data. The reference standard used during validation and study sample analysis should be obtained from an authentic and traceable source. The reference standard should be identical to the analyte. If this is not possible, an established form (e.g., salt or hydrate) of known quality should be used.

Suitable reference standards include compendial standards, commercially available standards or sufficiently characterized standards prepared in-house or by an external organization. A certificate of analysis (CoA) or an equivalent alternative is recommended to ensure quality and to provide information on the purity, storage conditions, retest/expiration date and batch number of the reference standard.

A CoA is not required for the IS as long as the suitability for use is demonstrated, e.g., a lack of analytical interference is shown for the substance itself or any impurities thereof.

When MS detection is used, the use of the stable isotope-labeled analyte as the IS is recommended whenever possible. However, it is important that the labeled standard is of high isotope purity and that no isotope exchange reaction occurs. The presence of unlabeled analyte should be checked and if unlabeled analyte is detected, the potential influence should be evaluated during method validation.

Stock and working solutions should only be prepared from reference standards that are within the stability period as documented in the CoA (either expiration date or the retest date).

B. Validation (3.2)

1. Selectivity (3.2.1)

Selectivity is the ability of an analytical method to differentiate and measure the analyte in the presence of potential interfering substances in the blank biological matrix.

Selectivity should be evaluated using blank samples (matrix samples processed without addition of an analyte or IS) obtained from at least six individual sources/lots (non-hemolyzed and non-lipemic). Use of fewer sources may be acceptable in the case of rare matrices. Selectivity for the IS should also be evaluated.

The evaluation of selectivity should demonstrate that no significant response attributable to interfering components is observed at the retention time(s) of the analyte or the IS in the blank samples. Responses attributable to interfering components should not be more than 20% of the

Contains Nonbinding Recommendations

analyte response at the LLOQ and not more than 5% of the IS response in the LLOQ sample for each matrix.

For the investigation of selectivity in lipemic matrices at least one source of matrix should be used. To be scientifically meaningful, the matrix used for these tests should be representative as much as possible of the expected study samples. A naturally lipemic matrix with abnormally high levels of triglycerides should be obtained from donors. Although it is recommended to use lipemic matrix from donors, if this is difficult to obtain, matrix can be spiked with triglycerides even though it may not be representative of study samples. However, if the drug impacts lipid metabolism or if the intended patient population is hyperlipidemic, the use of spiked samples is discouraged. This evaluation is not necessary for nonclinical studies unless the drug impacts lipid metabolism or is administered in a particular animal strain that is hyperlipidemic.

For the investigation of selectivity in hemolyzed matrices at least one source of matrix should be used. Hemolyzed matrices should be obtained by spiking matrix with hemolyzed whole blood (at least 2% V/V) to generate a visibly detectable hemolyzed sample.

2. *Specificity (3.2.2)*

Specificity is the ability of a bioanalytical method to detect and differentiate the analyte from other substances, including its related substances (e.g., substances that are structurally similar to the analyte, metabolites, isomers, impurities, degradation products formed during sample preparation, or concomitant medications that are expected to be used in the treatment of patients with the intended indication).

If the presence of related substances is anticipated in the biological matrix of interest, the impact of such substances should be evaluated during method validation, or alternatively, in the pre-dose study samples. In the case of LC-MS based methods, to assess the impact of such substances, the evaluation may include comparing the molecular weight of a potential interfering related substance with the analyte and chromatographic separation of the related substance from the analyte.

Responses detected and attributable to interfering components should not be more than 20% of the analyte response at the LLOQ and not more than 5% of the IS response in the LLOQ sample.

The possibility of back-conversion of a metabolite into the parent analyte during the successive steps of the analysis (including extraction procedures or in the MS source) should also be evaluated when relevant (e.g., potentially unstable metabolites such as ester analytes to ester/acidic metabolites, unstable N-oxides or glucuronide metabolites, lactone-ring structures). It is acknowledged that this evaluation will not be possible in the early stages of drug development of a new chemical entity when the metabolism is not yet evaluated. However, this issue should be investigated, and partial validation performed, if needed. The extent of back-

Contains Nonbinding Recommendations

conversion, if any, should be established and the impact on the study results should be discussed in the Bioanalytical Report.

3. *Matrix Effect (3.2.3)*

A matrix effect is defined as an alteration of the analyte response due to interfering and often unidentified component(s) in the sample matrix. During method validation the matrix effect between different independent sources/lots should be evaluated.

The matrix effect should be evaluated by analyzing at least three replicates of low and high QCs, each prepared using matrix from at least six different sources/lots. For each individual matrix sources/lots evaluated, the accuracy should be within $\pm 15\%$ of the nominal concentration and the precision (percent coefficient of variation (%CV)) should not be greater than 15%. Use of fewer sources/lots may be acceptable in the case of rare matrices.

The matrix effect should also be evaluated in relevant patient populations or special populations (e.g., hepatically impaired or renally impaired) when available. An additional evaluation of the matrix effect is recommended using hemolyzed or lipemic matrix samples during method validation on a case-by-case basis, especially when these conditions are expected to occur within the study.

4. *Calibration Curve and Range (3.2.4)*

The calibration curve demonstrates the relationship between the nominal analyte concentration and the response of the analytical platform to the analyte. Calibration standards, prepared by spiking matrix with a known quantity of analyte(s), span the calibration range and comprise the calibration curve. Calibration standards should be prepared in the same biological matrix as the study samples. The calibration range is defined by the LLOQ, which is the lowest calibration standard, and the ULOQ, which is the highest calibration standard. There should be one calibration curve for each analyte studied during method validation and for each analytical run.

A calibration curve should be generated with a blank sample, a zero sample (blank sample spiked with IS), and at least six concentration levels of calibration standards, including the LLOQ and the ULOQ.

A simple regression model that adequately describes the concentration-response relationship should be used. The selection of the regression model should be directed by written procedures. The regression model, weighting scheme, and transformation should be determined during the method validation. Blank and zero samples should not be included in the determination of the regression equation for the calibration curve. Each calibration standard may be analyzed in replicate, in which case data from all acceptable replicates should be used in the regression analysis.

The calibration curve parameters should be reported (e.g., slope and intercept in the case of a linear model). The back-calculated concentrations of the calibration standards should be presented together with the calculated mean accuracy and precision values. All acceptable curves

Contains Nonbinding Recommendations

obtained during validation, based on a minimum of three independent runs over several days, should be reported. The accuracy of the back-calculated concentrations of each calibration standard should be within $\pm 20\%$ of the nominal concentration at the LLOQ and within $\pm 15\%$ at all the other levels. At least 75% of the calibration standards with a minimum of six calibration standard levels should meet the above criteria.

In the case that replicates are used, the criteria (within $\pm 15\%$ or $\pm 20\%$ for LLOQ) should also be fulfilled for at least 50% of the calibration standards tested per concentration level. In the case that a calibration standard does not comply with these criteria, this calibration standard sample should be rejected, and the calibration curve without this calibration standard should be reevaluated, including regression analysis. For accuracy and precision runs, if all replicates of the LLOQ or the ULOQ calibration standard in a run are rejected, then the run should be rejected, the possible source of the failure should be determined and the method revised, if necessary. If the next validation run also fails, then the method should be revised before restarting validation.

The calibration curve should be prepared using freshly spiked calibration standards in at least one assessment. Subsequently, frozen calibration standards can be used within their defined period of stability.

5. *Accuracy and Precision (3.2.5)*

a. *Preparation of quality control samples (3.2.5.1)*

The QCs are intended to mimic study samples and should be prepared by spiking matrix with a known quantity of analyte, storing them under the conditions anticipated for study samples and analyzing them to assess the validity of the analytical method.

Calibration standards and the QCs should be prepared from separate stock solutions in order to avoid biased estimations which are not related to the analytical performance of the method. If calibration standards and the QCs may be prepared from the same stock solution, the accuracy and stability of the stock solution should be verified. A single source of blank matrix may be used, which should be free of interference or matrix effects, as described in section III.B.3 (3.2.3).

During method validation the QCs for accuracy and precision runs should be prepared at a minimum of 4 concentration levels within the calibration curve range: the LLOQ, within 3 times of the LLOQ (low QC), around 30% to 50% of the calibration curve range (medium QC) and at least 75% of the ULOQ (high QC).

For non-accuracy and precision validation runs, low, medium, and high QCs may be analyzed in duplicate. These QCs, along with the calibration standards, will provide the basis for the acceptance or rejection of the run.

b. *Evaluation of accuracy and precision (3.2.5.2)*

Contains Nonbinding Recommendations

Accuracy and precision should be determined by analyzing the QCs within each run (within-run) and in different runs (between-run). Accuracy and precision should be evaluated using the same runs and data.

Within-run accuracy and precision should be evaluated by analyzing at least five replicates at each QC concentration level in each analytical run. Between-run accuracy and precision should be evaluated by analyzing each QC concentration level in at least three analytical runs over at least two days. To enable the evaluation of any trends over time within one run, it is recommended to demonstrate accuracy and precision of the QCs over at least one of the runs in a size equivalent to a prospective analytical run of study samples. Reported method validation data and the determination of accuracy and precision should include all results obtained, including individual QCs outside of the acceptance criteria, except those cases where errors are obvious and documented. Within-run accuracy and precision data should be reported for each run. If the within-run accuracy or precision criteria are not met in all runs, an overall estimate of within-run accuracy and precision for each QC level should be calculated. Between-run (intermediate) precision and accuracy should be calculated by combining the data from all runs.

The calibration curves for these assessments should be prepared using freshly spiked calibration standards in at least one run. If freshly spiked calibration standards are not used in the other runs, stability of the frozen calibration standards should be demonstrated.

The accuracy at each concentration level should be within $\pm 15\%$ of the nominal concentration, except at the LLOQ, where it should be within $\pm 20\%$. The precision (%CV) of the concentrations determined at each level should not exceed 15%, except at the LLOQ, where it should not exceed 20%. For non-accuracy and precision validation runs, at least 2/3 of the total QCs and at least 50% at each concentration level should be within $\pm 15\%$ of the nominal values.

6. *Carryover (3.2.6)*

Carryover is an alteration of a measured concentration due to residual analyte from a preceding sample that remains in the analytical instrument.

Carryover should be assessed and minimized during method development. During validation carryover should be assessed by analyzing blank samples after the calibration standard at the ULOQ. Carryover in the blank samples following the highest calibration standard should not be greater than 20% of the analyte response at the LLOQ and 5% of the response for the IS. If it appears that carryover is unavoidable, study samples should not be randomized. Specific measures should be considered, validated, and applied during the analysis of the study samples, so that carryover does not affect accuracy and precision. This could include the injection of blank sample(s) after samples with an expected high concentration, before the next study sample.

7. *Dilution Integrity (3.2.7)*

Dilution integrity is the assessment of the sample dilution procedure, when required, to confirm that it does not impact the accuracy and precision of the measured concentration of the analyte.

Contains Nonbinding Recommendations

The same matrix from the same species used for preparation of the QCs should be used for dilution.

Dilution QCs should be prepared with analyte concentrations in matrix that are greater than the ULOQ and then diluted with blank matrix. At least five replicates per dilution factor should be tested in one run to determine if concentrations are accurately and precisely measured within the calibration range. The dilution factor(s) and concentrations applied during study sample analysis should be within the range of the dilution factors and concentrations evaluated during validation. The mean accuracy of the dilution QCs should be within $\pm 15\%$ of the nominal concentration and the precision (%CV) should not exceed 15%.

In the cases of rare matrices, use of a surrogate matrix for dilution may be acceptable. It should be demonstrated that this does not affect precision and accuracy.

8. *Stability (3.2.8)*

Stability evaluations should be carried out to ensure that every step taken during sample preparation, processing, and analysis as well as the storage conditions used do not affect the concentration of the analyte.

The storage and analytical conditions applied to the stability tests, such as the sample storage times and temperatures, sample matrix, anticoagulant and container materials, should reflect those used for the study samples. Reference to data published in the literature is not considered sufficient. Validation of storage periods should be performed on QCs that have been stored for a time that is equal to or longer than the study sample storage periods.

Stability of the analyte in the matrix is evaluated using low and high concentration QCs. Aliquots of the low and high QCs are analyzed at time zero and after the applied storage conditions that are to be evaluated. One bulk QC should be prepared at each concentration level. For each concentration tested, the bulk sample should be divided into a minimum of three aliquots that will be stored, stressed, and analyzed.

The QCs should be analyzed against a calibration curve, obtained from freshly spiked calibration standards in a run with its corresponding freshly spiked QCs or QCs for which stability has been proven. The mean concentration at each QC level should be within $\pm 15\%$ of the nominal concentration. If the concentrations of the study samples are consistently higher than the ULOQ of the calibration range, the concentration of the high QC should be adjusted to reflect these higher concentrations. It is recognized that this may not be possible in nonclinical studies due to solubility limitations.

For fixed dose combination products and specifically labeled drug regimens, the freeze-thaw, bench top, and long-term stability tests of an analyte in matrix should be conducted with the matrix spiked with all of the dosed compounds.

Contains Nonbinding Recommendations

The following stability tests should be evaluated:

(1) Stability of the Analyte in Matrix

Freeze-Thaw Stability in Matrix

To assess the impact of repeatedly removing samples from frozen storage, the stability of the analyte should be assessed after multiple cycles of freezing and thawing. Low and high QCs should be thawed and analyzed according to the same procedures as the study samples. QCs should be kept frozen for at least 12 hours between the thawing cycles. QCs for freeze-thaw stability should be assessed using freshly prepared calibration standards and QCs, or QCs for which stability has been proven. The number of freeze-thaw cycles validated should equal or exceed that of the freeze-thaw cycles undergone by the study samples, but a minimum of three cycles should be conducted.

Bench Top (Short-Term) Stability in Matrix

Bench top matrix stability experiments should be designed and conducted to cover the laboratory handling conditions for the study samples.

Low and high QCs should be thawed in the same manner as the study samples and kept on the bench top at the same temperature and for at least the same duration as the study samples.

The total time on the bench top should be concurrent; additive exposure to bench top conditions should not be used (i.e., time from each freeze-thaw evaluation should not be added up).

Long-Term Stability in Matrix

The long-term stability of the analyte in matrix stored in the freezer should be established. Low and high QCs should be stored in the freezer under the same storage conditions and at least for the same duration as the study samples.

For chemical drugs, the stability at one temperature (e.g., -20°C) can be extrapolated to lower temperatures (e.g., $-70^{\circ}/-80^{\circ}\text{C}$).

For biological drugs, a bracketing approach can be applied, e.g., in the case that the stability has been demonstrated at $-70^{\circ}/-80^{\circ}\text{C}$ and at -20°C , then it is not necessary to investigate the stability at temperatures in between those two points at which study samples will be stored.

(2) Stability of the Analyte in Processed Samples

The stability of processed samples, including the time until completion of analysis (in the autosampler/instrument), should be determined. For example:

Contains Nonbinding Recommendations

- Stability of the processed sample under the storage conditions to be used during the analysis of study samples (dry extract or in the injection phase)
- On-instrument/autosampler stability of the processed sample at injector or autosampler temperature

The total time that a processed sample is stored should be concurrent (i.e., autosampler and other storage times should not be added together).

(3) Stability of the Analyte and IS in Stock and Working Solutions

The stability of the stock and working solutions of the analyte and IS should be determined under the storage conditions used during the analysis of study samples by using the lowest and the highest concentrations of these solutions. They should be assessed using the response of the detector. Stability of the stock and working solutions should be tested with an appropriate dilution, taking into consideration the linearity and measuring range of the detector. If the stability varies with concentration, then the stability of all concentrations of the stock and working solutions should be assessed. If no isotopic exchange occurs for the stable isotopically labeled IS under the same storage conditions as the analyte for which the stability is demonstrated, then no additional stability determinations for the IS are necessary. If the reference standard expires, or it is past the retest date, the stability of the stock solutions made previously with this lot of reference standard are defined by the expiration or retest date established for the stock solution. Stock and working solutions should not be made from reference standards solely for extending the expiration date for the use of the reference standard.

In addition, the following test should be performed, if applicable:

(4) Stability of the Analyte in Whole Blood

Sufficient attention should be paid to the stability of the analyte in the sampled matrix (blood) directly after collection from subjects and prior to preparation for storage to ensure that the concentrations obtained by the analytical method reflect the concentrations of the analyte in the subject's blood at the time of sample collection.

If the matrix used is plasma, the stability of the analyte in blood should be evaluated during method development (e.g., using an exploratory method in blood) or during method validation. The results should be provided in the Validation Report.

9. *Reinjection Reproducibility (3.2.9)*

Reproducibility of the method is assessed by replicate measurements of the QCs and is usually included in the assessment of precision and accuracy. However, if samples could be reinjected (e.g., in the case of instrument interruptions or other reasons such as equipment failure), reinjection reproducibility should be evaluated to establish the viability of the processed samples and to support their storage prior to reinjection.

Contains Nonbinding Recommendations

Reinjection reproducibility is assessed by reinjecting a run that is comprised of calibration standards and a minimum of five replicates of the low, middle, and high QCs after storage. The precision and accuracy of the reinjected QCs establish the viability of the processed samples.

The results should be included in the Validation Report or provided in the Bioanalytical Report of the study where it was conducted.

C. Study Sample Analysis (3.3)

The analysis of study samples can be carried out after validation has been completed, however, it is understood that some parameters may be completed at a later stage (e.g., long-term stability). By the time the data are submitted to a regulatory authority, the bioanalytical method validation should have been completed. The study samples, QCs, and calibration standards should be processed in accordance with the validated analytical method. If system suitability is assessed, a predefined specific study plan, protocol, or SOP should be used. System suitability, including apparatus conditioning and instrument performance, should be determined using samples that are independent of the calibration standards and QCs for the run. Subject samples should not be used for system suitability. The IS responses of the study samples should be monitored to determine whether there is systemic IS variability. Refer to Table 1 for recommendations regarding documentation.

1. Analytical Run (3.3.1)

An analytical run should consist of a blank sample (processed matrix sample without analyte and without IS), a zero sample (processed matrix with IS), calibration standards at a minimum of six concentration levels, at least three levels of QCs (low, medium, and high) in duplicate (or at least 5% of the number of study samples, whichever is higher) and the study samples to be analyzed. The QCs should be interspersed in the run in such a way that the accuracy and precision of the whole run is ensured. Study samples should always be bracketed by QCs.

The calibration standards and QCs should be spiked independently using separately prepared stock solutions unless the accuracy and stability of the stock solutions have been verified. All samples (calibration standards, QCs, and study samples) should be processed and extracted as one single batch of samples in the order in which they are intended to be analyzed. Analyzing samples that were processed as several separate batches in a single analytical run is discouraged. If such an approach cannot be avoided, for instance due to bench top stability limitations, each batch of samples should include low, medium, and high QCs.

For comparative BA/BE studies, it is advisable to analyze all samples of one subject together in one analytical run to reduce variability.

The impact of any carryover that occurs during study sample analysis should be assessed and reported (refer to section III.B.6 (3.2.6)). If carryover is detected, its impact on the measured concentrations should be mitigated (e.g., non-randomization of study samples, injection of blank

Contains Nonbinding Recommendations

samples after samples with an expected high concentration) or the validity of the reported concentrations should be justified in the Bioanalytical Report.

2. Acceptance Criteria for an Analytical Run (3.3.2)

Criteria for the acceptance or rejection of an analytical run should be defined in the protocol, in the study plan, or in an SOP. In the case that a run contains multiple batches, acceptance criteria should be applied to the whole run and to the individual batches. It is possible for the run to meet acceptance criteria, even if a batch within that run is rejected for failing to meet the batch acceptance criteria. Calibration standards in a failed batch should not be used to support the acceptance of other batches within the analytical run.

The back-calculated concentrations of the calibration standards should be within $\pm 15\%$ of the nominal value, except for the LLOQ for which it should be within $\pm 20\%$. At least 75% of the calibration standard concentrations, which should include a minimum of six concentration levels, should fulfill these criteria. If more than six calibration standard levels are used and one of the calibration standards does not meet the criteria, this calibration standard should be rejected and the calibration curve without this calibration standard should be reevaluated, and a new regression analysis performed.

If the rejected calibration standard is the LLOQ, the new lower limit for this analytical run should be the next lowest acceptable calibration standard of the calibration curve. This new lower limit calibration standard should retain its original acceptance criteria (i.e., $\pm 15\%$). If the highest calibration standard is rejected, the ULOQ for this analytical run should be the next acceptable highest calibration standard of the calibration curve. The revised calibration range should cover at least three QC concentration levels (low, medium, and high). Study samples outside of the revised range should be reanalyzed. If replicate calibration standards are used and only one of the LLOQ or ULOQ standards fails, the calibration range is unchanged.

At least 2/3 of the total QCs and at least 50% at each concentration level should be within $\pm 15\%$ of the nominal values. If these criteria are not fulfilled the analytical run should be rejected. A new analytical batch should be prepared for all study samples within the failed analytical run for subsequent analysis. In the cases where the failure is due to an assignable technical cause, samples may be reinjected.

Analytical runs containing samples that are diluted and reanalyzed should include dilution QCs to verify the accuracy and precision of the dilution method during study sample analysis. The concentration of the dilution QCs should exceed that of the study samples being diluted (or of the ULOQ) and they should be diluted using the same dilution factor. If multiple dilution factors are used in one analytical run, then dilution QCs need only be diluted by the highest and lowest dilution factors. The within-run acceptance criteria of the dilution QC(s) will only affect the acceptance of the diluted study samples and not the outcome of the analytical run.

When several analytes are assayed simultaneously, there should be one calibration curve for each analyte studied. If an analytical run is acceptable for one analyte but has to be rejected for another analyte, the data for the accepted analyte should be used. The determination of the

Contains Nonbinding Recommendations

rejected analyte requires reextraction and analysis only for the analyte that is reanalyzed. Only data for this reanalyzed analyte needs to be reported.

The back-calculated concentrations of the calibration standards and QCs of passed and accepted runs should be reported. The overall (between-run) accuracy and precision of the QCs of all accepted runs should be calculated at each concentration level and reported in the analytical report (refer to section VIII. (8) Documentation and Table 1). If the overall mean accuracy and/or precision fails the 15% criterion, an investigation to determine the cause of the deviation should be conducted. In the case of comparative BA/BE studies, it may result in the rejection of the data.

3. Calibration Range (3.3.3)

If a narrow range of analyte concentrations of the study samples is known or anticipated before the start of study sample analysis, it is recommended to either narrow the calibration curve range, adapt the concentrations of the QCs, or add new QCs at different concentration levels as appropriate, to adequately reflect the concentrations of the study samples.

At the intended therapeutic dose(s), if an unanticipated clustering of study samples at one end of the calibration curve is encountered after the start of sample analysis, the analysis should be stopped and either the standard calibration range narrowed (i.e., partial validation), existing QC concentrations revised, or QCs at additional concentrations added to the original curve within the observed range before continuing with study sample analysis. It is not necessary to reanalyze samples analyzed before optimizing the calibration curve range or QC concentrations.

The same applies if a large number of the analyte concentrations of the study samples are above the ULOQ. The calibration curve range should be changed, if possible, and QC(s) added, or their concentrations modified. If it is not possible to change the calibration curve range or the number of samples with a concentration above the ULOQ is not large, samples should be diluted according to the validated dilution method.

At least two QC levels should fall within the range of concentrations measured in study samples. If the calibration curve range is changed, the bioanalytical method should be revalidated (partial validation) to verify the response function and to ensure accuracy and precision.

4. Reanalysis of Study Samples (3.3.4)

Possible reasons for reanalysis of study samples, the number of replicates, and the decision criteria to select the value to be reported should be predefined in the protocol, study plan, or SOP before the actual start of the analysis of the study samples. For study samples in which multiple analytes are being analyzed, a valid result for one analyte should not be rejected if the other analyte fails the acceptance criteria.

The number of samples (and percentage of total number of samples) that have been reanalyzed should be reported and discussed in the Bioanalytical Report. For comparative BA/BE studies, a separate table should report values from rejected runs.

Contains Nonbinding Recommendations

Some examples of reasons for study sample reanalysis are:

- Rejection of an analytical run because the run failed the acceptance criteria with regard to accuracy of the calibration standards and/or the precision and accuracy of the QCs
- IS response significantly different from the response for the calibration standards and QCs (as predefined in an SOP)
- The concentration obtained is above the ULOQ
- The concentration observed is below the revised LLOQ in runs where the lowest calibration standard has been rejected from a calibration curve, resulting in a higher LLOQ compared with other runs
- Improper sample injection or malfunction of equipment
- The diluted study sample is below the LLOQ
- Identification of quantifiable analyte levels in pre-dose samples, control, or placebo samples
- Poor chromatography (as predefined in an SOP)

For comparative BA/BE studies, reanalysis of study samples for a PK reason (e.g., a sample concentration does not fit with the expected profile) should not be done, as it may bias the study result.

Any reanalyzed samples should be identified in the Bioanalytical Report and the initial value, the reason for reanalysis, the values obtained in the reanalyzes, the final accepted value and a justification for the acceptance should be provided. Further, a summary table of the total number of samples that have been reanalyzed for each reason should be provided. In cases where the first analysis yields a nonreportable result, a single reanalysis is considered sufficient (e.g., concentration above the ULOQ or equipment malfunction). In cases where the value needs to be confirmed (e.g., pre-dose sample with measurable concentrations) replicate determinations are required if sample volume allows.

The safety of trial subjects should take precedence over any other aspect of the trial. Consequently, there may be other circumstances when it is necessary to reanalyze specific study samples for the purpose of a safety investigation.

5. Reinjection of Study Samples (3.3.5)

Reinjection of processed samples can be made in the case of equipment failure if reinjection reproducibility has been demonstrated during validation or provided in the Bioanalytical Report where it was conducted. Reinjection of a full analytical run or of individual calibration standards

Contains Nonbinding Recommendations

or QCs simply because the calibration standards or QCs failed, without any identified analytical cause, should not be done.

6. *Integration of Chromatograms (3.3.6)*

Chromatogram integration and reintegration should be described in a study plan, protocol, or SOP. Any deviation from the procedures described *a priori* should be discussed in the Bioanalytical Report. The list of chromatograms that required reintegration, including any manual integrations, and the reasons for reintegration should be included in the Bioanalytical Report. Original and reintegrated chromatograms and initial and repeat integration results should be kept for future reference and submitted in the Bioanalytical Report for comparative BA/BE studies.

IV. LIGAND BINDING ASSAYS (4)

A. Key Reagents (4.1)

1. *Reference Standard (4.1.1)*

The reference standard should be well characterized and documented (e.g., CoA and origin). A biological drug has a highly complex structure and its reactivity with binding reagents for bioanalysis may be influenced by a change in the manufacturing process of the drug substance. It is recommended that the manufacturing batch of the reference standard used for the preparation of calibration standards and QCs is derived from the same batch of drug substance as that used for dosing in the nonclinical and clinical studies whenever possible. If the reference standard batch used for bioanalysis is changed, bioanalytical evaluation should be carried out with QCs from the original material and the new material prior to use to ensure that the performance characteristics of the method are within the acceptance criteria.

2. *Critical Reagents (4.1.2)*

Critical reagents, including binding reagents (e.g., binding proteins, aptamers, antibodies, or conjugated antibodies) and those containing enzymatic moieties, have direct impact on the results of the assay and, therefore, their quality should be assured. Critical reagents bind the analyte and, upon interaction, lead to an instrument signal corresponding to the analyte concentration. The critical reagents should be identified and defined in the assay method.

Reliable procurement of critical reagents, whether manufactured in-house or purchased commercially, should be considered early in method development. The data sheet for the critical reagent should include at a minimum identity, source, batch/lot number, purity (if applicable), concentration (if applicable), and stability/retest date/storage conditions (refer to Table 1). Additional characteristics may be warranted.

A critical reagent lifecycle management procedure is necessary to ensure consistency between the original and new batches of critical reagents. Reagent performance should be evaluated using

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的注册商标或商标。本书中有未标明 TM 标志和 © 标志之处。