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Group Condition: [Quantitative][PM-Test-STD]

File Help

Group Condition

- Quantitative Information
 - Element Information
 - Measurement Condition
 - Calibration Curve
 - Sensitivity Coefficient
 - BG/Intensity Calculati
 - Overlap Correction
 - Internal Standard
 - Result Format
 - L/H Decision
 - Additional Calculation
 - Sensitivity Coefficient
 - Meas. Condition List
- Standard Sample
 - Standard Sample
 - Standardization
 - Drift Check
 - Master Curve
 - Blank Sample

Group Measurement Condition

Aperture: 30 Atmosphere: Vac.

Spin Slow Evac./Air Pu

Rotation: 0 Vacuum Stabiliz

Pressure to: 30 Pa

Optimization

Online Output

Quant: Calculat

Output: Print File

Display: Cond. Intern. Each Quant.

Matching (No Value)

Information of a Sample

Sample Form: Bulk Film Filter Pap

Compound Form: Metal Oxide Reference System: default

Elements Included

Layer Structure:

1	Layer1	Pb	Ba	Sr	Br	Se	As	Sb	Sn	Cd	Zn
		Cu	Ni	Co	Fe	Mn	Cr	V	Ti	So	Ca
		K	Cl	S	P	Si	Al	Mg	Na		

Group Condition

- Quantitative Information
 - Element Information
 - Measurement Condition
 - Calibration Curve
 - Sensitivity Coefficient
 - BG/Intensity Calculati
 - Overlap Correction
 - Internal Standard
 - Result Format
 - L/H Decision
 - Additional Calculation
 - Sensitivity Coefficient
 - Meas. Condition List
- Standard Sample
 - Standard Sample
 - Standardization
 - Drift Check
 - Master Curve
 - Blank Sample

Layer	Info	Analyte	Fixed.Val.	Unit	Proc.-Calc.	Line	Ms	BG	Ds	Bo	Ov
1	Layer1			mg/cm2	Total		x	x	x	x	x
1	Elem. Pb			ug/cm2	Quant-EC PbLb1		o	x	x	x	x
1	Elem. Ba			ug/cm2	Quant-EC BaLa		o	x	x	x	x
1	Elem. Sr			ug/cm2	Quant-EC SrKa		o	x	x	x	x
1	Elem. Br			ug/cm2	Quant-EC BrKa		o	x	x	x	x
1	Elem. Se			ug/cm2	Quant-EC SeKa		o	x	x	x	x
1	Elem. As			ug/cm2	Quant-EC AsKa		o	x	x	x	x

Element Information

[Layer 1] Layer Info.

Compound Name:

Unit:

Proc.-Calc.

Quant. Correct Fix Flux Balance Total

Type of Quantification

FP BG-FP Calibration Cu

Fixed Value

Fix: 0.0000

Input: Fixed

Bg Analysis:

Prev Next

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Group Condition

- Quantitative Information
 - Element Information
 - Measurement Condition
 - Calibration Curve
 - Sensitivity Coefficient
 - BG/Intensity Calculati
 - Overlap Correction
 - Internal Standard
 - Result Format
 - L/H Decision
 - Additional Calculation
 - Sensitivity Coefficient
 - Meas. Condition List
- Standard Sample
 - Standard Sample
 - Standardization
 - Drift Check
 - Master Curve
 - Blank Sample

Layer	Info	Analyte	Fixed.Val.	Unit	Proc.-Calc.	Line	Ms	BG	Ds	Bo	Ov
1	Layer1			mg/cm2	Total		x	x	x	x	x
1	Elem. Pb			ug/cm2	Quant-EC PbLb1		o	x	x	x	x
1	Elem. Ba			ug/cm2	Quant-EC BaLa		o	x	x	x	x
1	Elem. Sr			ug/cm2	Quant-EC SrKa		o	x	x	x	x
1	Elem. Br			ug/cm2	Quant-EC BrKa		o	x	x	x	x
1	Elem. Se			ug/cm2	Quant-EC SeKa		o	x	x	x	x
1	Elem. As			ug/cm2	Quant-EC AsKa		o	x	x	x	x

Element Information

[Layer 1] Elem. Info.

Compound Name:

Unit:

Proc.-Calc.

Quant. Correct Fix Flux Balance Total

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
Group Condition: [Quantitative][PM-Test-STD]

File Help

Group Condition

- Quantitative Information
- Element Information
- Measurement Condition
- Calibration Curve
- Sensitivity Coefficient
- BG/Intensity Calculati
- Overlap Correction
- Internal Standard
- Result Format
- L/H Decision
- Additional Calculation
- Sensitivity Coefficient
- Meas. Condition List
- Standard Sample
- Standard Sample
- Standardization
- Drift Check
- Master Curve
- Blank Sample

Layer	Info	Analyte	Fixed.Val.	Unit	Proc.-Calc.	Line	Ms	BG	Ds	Be	Ov	Is
1 Layer1												
1 Layer												
1 Elem. Pb				mg/cm2	Total		x	x	x	x	x	x
1 Elem. Ba				ug/cm2	Quant.-EC	PbLb1	o	o	x	x	x	o
1 Elem. Sr				ug/cm2	Quant.-EC	SrKa	o	o	x	x	x	o
1 Elem. Br				ug/cm2	Quant.-EC	BrKa	o	o	x	x	x	o
1 Elem. Se				ug/cm2	Quant.-EC	SeKa	o	o	x	x	x	o

Measurement Condition

Do Measure

X-ray Generator: Target: Pb, Voltage (kV): 40, Current (mA): 70, Filter: None

Goniometer: Mode: Fix, Angle(deg): 28.200, Meas. Time(sec): 10.00, BG Correction: 1 point, Func.: No. Angle

No.	Angle
1	28.700 7 10.00

Total Meas. 20.0


Spectrometer: Slit: Standard, Attenuator: Off, Crystal: 1:LiF

Detector System: Detector: SC, PHA Low Level: 20, PHA High Level: 80

Auto Sensitivity Control: Control: mA, Attenuation: 1/20, Threshold(kcps): 300

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The image displays two screenshots from the Shimadzu software interface, illustrating the process of confirming the analysis angle.


Top Screenshot: Measurement Condition Dialog

- Measurement:** Position 1, Measurement 2theta. *1 选择进行 2θ 扫描*
- 2theta Condition:** Scanning Fixed Angle, Step Continuous, Peak Angle 62.97, Start 61, End Angle 65. *2 选择少步或连续*
- 3 设定范围包含上面峰的角度** (Set range to include the peak angle above)
- 4 选择步长角度 0.05 或 0.02** (Select step angle 0.05 or 0.02)

Bottom Screenshot: Profile Monitor and Measurement Dialog

- Profile Monitor:** Shows a graph of intensity [kcps] vs. angle [deg]. A peak is visible at approximately 62.5 degrees. A context menu is open over the peak with options: None, Set Range, and Set Ang. *Set Ang 如果峰值位置不对, 鼠标左右移动, 确定最佳位置*
- Measurement Dialog:** Position 8, Measurement 2theta, PHD, Intensit. *1 选择标样位置*
- Start Button:** Circled in red, indicating the start of the measurement.
- 2theta Condition:** Scanning Fixed Angle, Step Continuous, Peak Angle 62.97, Step Angle 0.020, Start 61, Speed 8.0, End Angle 65, Time 0.15.
- PHD Condition:** Peak Angle 62.97, Start 0, Step Width 2, End 200, Time 1.

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Group Condition: [Quantitative][PM-Test-STD]

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

Group Condition: [Quantitative][PM-Test-STD]

File Help

Group Condition

- Quantitative Information
 - Element Information
 - Measurement Condition
 - Calibration Curve
 - Sensitivity Coefficient
 - BG/Intensity Calculati
 - Overlap Correction
 - Internal Standard
 - Result Format
 - L/H Decision
 - Additional Calculation
 - Sensitivity Coefficient
- Standard Sample
 - Standard Sample
 - Standardization
 - Drift Check
 - Master Curve
 - Blank Sample

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Select Group

[Analysis] : [Ready]

File Edit View Test Measurement Help

Sample Name... Process Data...

Schedule of Analysis

No	Sample Name	Purpose of Measurement
1	[Instrument Setup]Shutdown	
1	Keeps Xray	
2	[Quantitative]PM-test1	-F-
1	Blank	Std.--Intensity
2	33964	Std.--Intensity
3	33951	Std.--Intensity
4	33953	Std.--Intensity
5	33952	Std.--Intensity
6	33946	Std.--Intensity
7	33949	Std.--Intensity
8	33950	Std.--Intensity

Sample Registration

Analysis Instrument Setup PH

Analytical Group: Easy Analysis...

[Quantitative]PM-test1

Detail Settings

Repeat: Pause:

Sample

Sample Name	Position
Blank	1
33964	2
33951	3
33953	4
33952	5
33946	6
33949	7
33950	8
33961	1
33948	2

Select Analytical Group


Easy Analysis Qual-Quant. Quantitative

Group Name	Date	Comment
薄膜吸附量-过零点-Pb	2020/11/11 11:11:51	Comment
薄膜吸附量-3mm	2019/03/15 17:15:22	Comment
薄膜吸附量-新-过零点	2019/03/15 12:15:28	Comment
薄膜吸附量-新1	2017/12/18 15:35:44	Comment
薄膜吸附量-guangxi2	2017/06/23 10:03:36	special for NaCl
薄膜吸附量-CaPbTi	2016/05/17 16:13:44	Comment
薄膜吸附量-guangxi1	2016/05/16 16:03:30	Comment
normal el@ant-SY	2016/04/18 16:06:08	Comment
PM-test1	2016/04/18 10:19:30	Comment
薄膜吸附量-加膜	2016/03/30 08:03:02	Comment
薄膜吸附量	2015/12/22 10:16:34	Comment

Purpose of Measurement

Standard Intensity Registration Coefficient Renewal

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Group Condition: [Quantitative][PM-Test-STD]

Layer	Info	Analyte	Fixed Val.	Unit	Proc.	Calc.	Line	Ms	BG	Ds	Bo
1	Layer1			Total	x	x	x	x	x	x	x
1	Elem.	Pb		ug/cm2	Quant.-EC	PbLb1	o	o	x	x	x
1	Elem.	Ba		ug/cm2	Quant.-EC	BaLa	o	o	x	x	x
1	Elem.	Sr		ug/cm2	Quant.-EC	SrKa	o	o	x	x	x
1	Elem.	Br		ug/cm2	Quant.-EC	BrKa	o	o	x	x	x
1	Elem.	Se		ug/cm2	Quant.-EC	SeKa	o	o	x	x	x
1	Elem.	As		ug/cm2	Quant.-EC	AsKa	o	o	x	x	x

Calibration Curve

Graph 1: X : Std Value Y : Measured In

Curve No. 1 Multiple Calibration Curve

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
Sample Registration

Analysis Instrument Setup PM Adjustment

Analytical Group: Easy Analysis...

[Quantitative]PM-test1 Unknown Output -F-

Repeat: ___ Pause: ___ New Page: ___ Detail Settings...

Sample Name	Position	Comment	Fixed Value	Special Analysis
1: 123	1	Comment	Unknown	
2: 124	2	Comment	Unknown	
3: 125	3	Comment	Unknown	
4:	4	Comment	Unknown	
5:	5	Comment	Unknown	
6:	6	Comment	Unknown	
7:	7	Comment	Unknown	
8:	8	Comment	Unknown	
9:	1	Comment	Unknown	
10:	2	Comment	Unknown	

Clear Apply Close

Select Analytical Group

Easy Analysis Qual-Quant. Quantitative

Group Name	Date	Comment
薄层吸附-过零点-Pb	2020/11/11 11:11:51	Comment
薄层吸附-3m	2019/03/15 17:15:22	Comment
薄层吸附-新-过零点	2019/03/15 12:15:28	Comment
薄层吸附-新1	2017/12/18 15:35:44	Comment
薄层吸附-新-guangpi2	2017/06/23 10:03:36	special for NaCl
薄层吸附-新-CaPb1	2016/05/17 16:19:44	Comment
薄层吸附-新-guangpi1	2016/05/16 16:03:30	Comment
Normal element-SY	2016/04/18 16:08:08	Comment
薄层吸附-新-加膜	2016/03/30 09:03:02	Comment
薄层吸附-新	2015/12/22 10:16:34	Comment

Purpose of Measurement: Unknown Intensity Coefficient

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