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- [1] Agency for Toxic Substances and Disease Registry, "Per- and Polyfluoroalkyl Substances (PFAS) and Your Health," 31 October 2018. [Online]. Available: <https://www.atsdr.cdc.gov/pfas/>. [Accessed 11 December 2018].
- [2] C. Lau, K. Anitole, C. Hodes, D. Lai, A. Pfahles-Hutchens and J. Seed, "Perfluoroalkyl Acids: A Review of Monitoring and Toxicological Findings," *Toxicological Sciences*, vol. 99, no. 2, pp. 366-394, 1 October 2007.
- [3] Environmental Protection Agency, "Lifetime Health Advisories and Health Effects Support Documents for Perfluorooctanoic Acid and Perfluorooctane Sulfonate," *Federal Register*, vol. 81, no. 101, pp. 33250-33251, 2016.
- [4] United States Environmental Protection Agency (US EPA), "PFOA, PFOS and Other PFASs: Basic Information about PFAS," 20 August 2018. [Online]. Available: <https://www.epa.gov/pfas/basic-information-pfas#exposed>. [Accessed 5 December 2018].
- [5] Sara Sahlin, Swedish Environmental Protection Agency, "PFAS in the Baltic Sea Region: Inventory of Awareness, Actions and Strategies Related to Highly Fluorinated Substances PFAS, including PFOS," 2017.
- [6] W. H. O. (WHO), Keeping our water clean: the case of water contamination in the Veneto Region, Italy, Copenhagen, Denmark: World Health Organization (WHO) Regional Office for Europe, 2016.
- [7] K. M. Rappazzo, E. Coffman and E. P. Hines, "Exposure to Perfluorinated Alkyl Substances and Health Outcomes in Children: A Systematic Review of the Epidemiologic Literature," *International Journal of Environmental Research and Public Health*, vol. 14, no. 7, p. 691, 2017.
- [8] United States Environmental Protection Agency, "Supporting Documents for Drinking Water Health Advisories for PFOA and PFOS," November 2016. [Online]. Available: <https://www.epa.gov/ground-water-and-drinking-water/supporting-documents-drinking-water-health-advisories-pfoa-and-pfos>. [Accessed 4 December 2018].
- [9] United States Environmental Protection Agency (US EPA), "U.S. State Resources about PFAS," [Online]. Available: <https://www.epa.gov/pfas/us-state-resources-about-pfas>. [Accessed 4 December 2018].
- [10] C. Alder, "Analysis of State-by-State Differences in PFAS Regulation," Northeastern University Social Science Environmental Health Research Institute, 2 October 2018. [Online]. Available: <https://pfasproject.com/2018/10/02/analysis-of-state-by-state-differences-in-pfas-regulation/>. [Accessed 11 December 2018].
- [11] California Water Boards, State Water Resources Control Board, "Perfluorooctanoic acid (PFOA) and Perfluorooctanesulfonic acid (PFOS)," 13 July 2018. [Online]. Available: https://www.waterboards.ca.gov/drinking_water/certlic/drinkin_gwater/PFOA_PFOS.html. [Accessed 13 December 2018].
- [12] X. C. Hu, D. Q. Andrews, A. B. Lindstrom, T. A. Bruton, L. A. Schaider, P. Grandjean, R. Lohmann, C. C. Carignan, A. Blum, S. A. Balan, C. P. Higgins and E. M. Sunderland, "Detection of Poly- and Perfluoroalkyl Substances (PFASs) in U.S. Drinking Water Linked to Industrial Sites, Military Fire Training Areas and Wastewater Treatment Plants," *Environmental Science & Technology Letters*, vol. 3, no. 10, pp. 344-350, 2016.

- [13] H. A. Kabore, S. Vo Duy, G. Munoz, M. Ladjji, M. Desrosiers, J. Liu, T. K. Sory and S. Sauve, "Worldwide drinking water occurrence and levels of newly-identified perfluoroalkyl and polyfluoroalkyl substances," *Science of The Total Environment*, Vols. 616-617, pp. 1089-1100, 2018.
- [14] U.S. Environmental Protection Agency, "US EPA Method 537: Determination of Selected Perfluorinated Alkyl Acids in Drinking Water by Solid Phase Extraction and Liquid Chromatography / Tandem Mass Spectrometry (LC/MS/MS)," Washington D.C., 2009.
- [15] United States Environmental Protection Agency, "US EPA - PFAS Research and Development," 14 August 2018. [Online]. Available: https://www.epa.gov/sites/production/files/2018-08/documents/r4_combined_presentations_.pdf. [Accessed 5 December 2018].
- [16] ASTM International, "ASTM D7979-17: Standard Test Method for Determination of Perfluorinated Compounds in Water, Sludge, Influent, Effluent and Wastewater by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS)," West Conshohocken, 2017.
- [17] ASTM International, "ASTM D7968-17a: Standard Test Method for Determination of Perfluorinated Compounds in Soil by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS)," West Conshohocken, 2017.
- [18] United States Environmental Protection Agency (US EPA), "Assessing and Managing Chemicals under TSCA, Fact Sheet: 2010/2015 PFOA Stewardship Program," [Online]. Available: <https://www.epa.gov/assessing-and-managing-chemicals-under-tsca/fact-sheet-20102015-pfoa-stewardship-program>. [Accessed 23 December 2018].
- [19] Secretariat of the Basel, Rotterdam and Stockholm Conventions (UNEP), "Two More Toxic Chemicals Recommended to be Eliminated as UN Scientific Committee Paves Way for Ban on Widely-Used PFOA [Press Release]," 18 September 2018. [Online]. Available: <http://www.brsmeas.org/Implementation/MediaResources/PressReleases/POPRC14PressReleases/tabid/7685/language/en-US/Default.aspx>. [Accessed 3 December 2018].
- [20] S. A. Meiburg, "Emerging Contaminants and Environmental Health," *North Carolina Medical Journal*, vol. 79, no. 5, pp. 315-316, 2018.
- [21] U.S. Environmental Protection Agency, "US EPA Method 537.1: Determination of Selected Per- and Polyfluorinated Alkyl Substances in Drinking Water by Solid Phase Extraction and Liquid Chromatography / Tandem Mass Spectrometry (LC/MS/MS)," Washington D.C., 2018.

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