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Fig.1 the HPLC-GC Chromatograms of MOSH and MOAH in the milk powder sample (RVEP 200612) from DRRR Proficiency Testing

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- [1] Foodwatch. International test of various canned baby milk products for their content of mineral oil hydrocarbons (MOSH/MOAH) [EB/OL]. (2019-10-24) [2020-02-04]. https://www.foodwatch.org/fileadmin/-INT/mineral_oil/documents/2019-10-24_Projectreport_babymilk_FINAL.pdf
- [2] European Food Safety Authority (EFSA), ARCELLA D, BAERT K, et al. Rapid risk assessment on the possible risk for public health due to the contamination of infant formula and follow-on formula by mineral oil aromatic hydrocarbons (MOAH)[J]. EFSA Supporting Publication. 2019, 16(11): -.
- [3] ARES. 2019. Summary of the roundtable workshop on the determination of MOAH in infant formula[C]. European Commission Joint Research Centre. 7564336 - 09/12/2019.
- [4] Food Federation Germany, Consumer Protection Consortium of the Federal States. Benchmark levels for mineral oil hydrocarbons (MOH) in foods [EB/OL]. (2022-10-12) [2022-11-10]. <https://www.lebensmittelverband.de/de/aktuell/20221012-aktualisierung-moh-orientierungswerte>.
- [5] Standing Committee on Plants, Animals, Food and Feed (SCoPAFF). Section Novel Food and Toxicological Safety of the Food Chain[EB/OL]. [2022-11-10]. https://ec.europa.eu/food/system/files/2022-05/reg-com_toxic_20220421_sum.pdf.
- [6] European Food Safety Authority (EFSA). Scientific opinion on mineral oil hydrocarbons in food[J]. EFSA Journal, 2012, 10(6): 2704.
- [7] 黄华, 武彦文, 李冰宁, 等. 食用植物油中矿物油污染物的分析进展 [J]. 中国油脂, 2018, 43(07): 97-101. [HUANG H, WU Y W, LI B N, et al. Progress in analysis of mineral oil contaminants in edible vegetable oils[J]. China Oils and Fats. 2018, 43(07): 97-101.
- [8] BRÜHL L. Occurrence, determination, and assessment of mineral oils in oilseeds and vegetable oils[J]. European Journal of Lipid Science and Technology[J]. 2016, 118(3): 361-372.
- [9] LIU L L, HUANG H, WU Y W, et al. Offline solid-phase extraction large volume injection gas chromatography for the analysis of mineral oil saturated hydrocarbons in commercial vegetable oils[J]. Journal of Oleo Science, 2017, 66(9): 981-990.

- [10] RUIZ J L H, LIÉBANAS J A, VIDAL J L M, et al. Offline solid-phase extraction and separation of mineral oil saturated hydrocarbons and mineral oil aromatic hydrocarbons in edible oils, and analysis via GC with a flame ionization detector[J]. *Foods*, 2021, 10(9): 2026.
- [11] GIÓMEZ-COCA R B, CERT R, PÉREZ-CAMINO M C, et al. Determination of saturated aliphatic hydrocarbons in vegetable oils[J]. *Grasasy Aceites*, 2016, 67(2): e127.
- [12] 李冰宁, 刘玲玲, 武彦文. 固相萃取 - 程序升温 - 气相色谱法与大容量固相萃取 - 气相色谱法测定食用油脂中饱和烃矿物油的比较 [J]. *分析实验室*, 2018, 37(6): 701-705. [LI B N, LIU L L, WU Y W. Comparison of solid phase extraction-programmed temperature vaporizer gas-chromatography and large scale-solid phase extraction-gas chromatography for the determination of mineral oil saturated hydrocarbons in edible oils and fats[J]. *Chinese Journal of Analysis Laboratory*, 2018, 37(6): 701-705.]
- [13] BARP L, KORNAUTH C, WUERGER T, et al. Mineral oil in human tissue, Part I: Concentrations and molecular mass distributions[J]. *Food and Chemical Toxicology*, 2014, 72: 312-321.
- [14] BIEDERMANN M, BARP L, KORNAUTH C, et al. Mineral oil in human tissues, Part II: Characterization of the accumulated hydrocarbons by comprehensive two-dimensional gas chromatography[J]. *Science of The Total Environment*, 2015, 506-507: 644-655.
- [15] BARP L, BIEDERMANN M, GROB K, et al. Accumulation of mineral oil saturated hydrocarbons (MOSH) in female Fischer 344 rats: Comparison with human data and consequences for risk assessment[J]. *Science of The Total Environment*, 2017, 275:1263-1278.
- [16] BIEDERMANN M, GROB K. On-line coupled high performance liquid chromatography-gas chromatography for the analysis of contamination by mineral oil. Part 1: Method of analysis[J]. *Journal of Chromatography A*, 2012, 1255: 56-75.
- [17] 武彦文, 李冰宁, 刘玲玲, 等. 基于高效液相色谱 - 气相色谱联用技术的矿物油分析研究进展 [J]. *分析化学*, 2021, 49(03): 341-349. [WU Y W, LI B N, LIU L L, et al. Research progress of mineral oil analysis based on high performance liquid chromatography-gas chromatography[J]. *Chinese Journal of Analytical Chemistry*, 2021, 49(03): 341-349.]
- [18] European Committee for Standardization (CEN) European standard (EN) 16995-2017 Foodstuffs-Vegetable oils and food-stuff on basis of vegetable oils-Determination of mineral oil saturated hydrocarbons (MOSH) and mineral oil aromatic hydrocarbons (MOAH) with on-line HPLC-GC-FID analysis[S]. Brussels: BSI standards Limited, 2017.
- [19] BRATINOVA S, HOEKSTRA E. Guidance on sampling, analysis and data reporting for the monitoring of mineral oil hydrocarbons in food and food contact materials[J]. Luxembourg: Publications Office of the European Union, 2019, ISBN 978-92-76-00172-0, doi:10.2760/208879, JRC 115694.
- [20] LIU L L, LI B N, OUYANG J, Wu Y W. Mild mixed-solvent extraction for determination of total mineral oil hydrocarbon contaminants in milk powder products [J]. *Food Chemistry*, 2020, 333: 127488.
- [21] 杨梦奇, 刘玲玲, 李冰宁, 等. 植物油基鱼罐头中矿物油污染物高灵敏检测方法的建立 [J]. *食品工业科技*, 2022, 43(12): 317-324. [YANG M Q, LIU L L, LI B N, et al. Establishment of Highly Sensitive Determination Method for Mineral Oil Contaminants in Vegetable Oil-based Canned Fish Products [J]. *Science and Technology of Food Industry*, 2022, 43(12): 317-324.]

- [22] German Society for Fat Science (DGF), Draft version for updating EN-16995:2017: October 2022-V4 Vegetable oils - Determination of mineral oil saturated hydrocarbons (MOSH) and aromatic hydrocarbons (MOAH) with online coupled HPLC-GC-FID analysis -Method for low limit of quantification.
- [23] 李鑫. 婴幼儿配方乳粉生产工艺、包装的现状研究 [J]. 中国乳业, 2021, 05: 85-90. [Xin L. Current Situation and Research on Production conditions and Packages[J]. China Dairy, 2021, 05: 85-90.]
- [24] 揭良, 苏米亚, 贾宏信, 等. 婴幼儿配方乳粉脂质母乳化研究进展 [J]. 乳业科学与技术, 2020, 43(3): 45-49. [JIE L, SU M Y, JIA H X, et al. Recent progress in infant formulas simulating the fat composition of human milk[J]. Journal of Dairy Science and Technology, 2020, 43(3): 45-49.]
- [25] BIEDERMANN M, FISELIER K, GROB K. Aromatic hydrocarbons of mineral oil origin in foods: method for determining the total concentration and first results[J]. Journal of Agricultural and Food chemistry, 2009, 57(19): 8711-8721.
- [26] NESTOLA M, SCHMIDT TC. Determination of mineral oil aromatic hydrocarbons in edible oils and fats by online liquid chromatography-gas chromatography-flame ionization detection - Evaluation of automated removal strategies for biogenic olefins[J]. Journal of Chromatography A, 2017, 1505: 69-76.
- [27] European Committee for Standardization (CEN) European standard (EN) 17517-2021 Animal feeding stuffs: Methods of sampling and analysis - Determination of mineral oil saturated hydrocarbons (MOSH) and mineral oil aromatic hydrocarbons (MOAH) with on-line HPLC-GC-FID analysis [S]. Brussels: BSI standards Limited, 2021.
- [28] Deutsches Institut für Normung, e.V. (DIN) SPEC 5010:2018-05 Testing of paper and board-Determination of the transfer of mineral oil hydrocarbons from food contact materials manufactured with portions of recycled pulp[S]. Berlin, 2018.
- [29] 中国合格评定国家认可委员会. CNAS-GL002-2018 能力验证结果的统计处理和评价指南 [S]. 北京: 中国合格评定国家认可委员会, 2018. [CNAS-GL002-2018 Guidance on Statistics Treatment of Proficiency Testing Results and Performance Evaluation[S]. Beijing: China National Accreditation Service for Conformity Assessment, 2018.]
- [30] 全国认证认可标准化技术委员会. GB/T 27043-2012 合格评定 能力验证的通用要求 [S]. 北京: 中国标准出版社, 2012. [GB 27043-2012 Conformity assessment—General requirements for proficiency testing[S]. Beijing: China Standard Press, 2012.]

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